

MONITEX

ROTAMIX

GX330

Alginate Mixer

Smooth mixture in seconds for a quick and accurate impression

Reduce wasted alginate powder

Easy to use and clean

Save space and time



MONITEX

ROTAMIX GX330 Alginate Mixer

The automated RotaMIX alginate mixer makes homogeneous and smoothly mixed alginate pastes to create the best impressions. With an alginate mixer in the clinic, dentists can be sure that there are no compromises, such as those arising from manual mixing, so that all procedures may be standardized.



Simple steps to creating a perfect impression:



1. Add alginate powder and water into the mixing cup.

2. Moisten with the stirrer for a few seconds and close the mixing cup.





3. Put the mixing cup into the mixing well.

4. Cover the lid and press "START".





5. After mixing, use the spatula and remove the mixed impression.

3 Default Mixing Modes

Scoops	Mixing Time
☪	8 secs
☪ ☪	10 secs
☪ ☪ ☪	12 secs

Specifications:

- Timer Range 1 ~ 16 seconds
- Mixing Modes 3 (Default: 8, 10, 12 secs)
- Speed 3,400 rpm
- Power Input AC 110 or 220V, 50/60 Hz
- Dimensions 250(L) x 205(W) x 290(H) mm
- Net Weight 17.5 kg
- Accessories Spatula/Stirrer/Storage Case/
Water Dosing Bottle/Mixing Cups(x4)



MONITEX INDUSTRIAL CO., LTD

5F, 70, Sec. 1, Guangfu Rd., Sanchong Dist., New Taipei City 24158, TAIWAN
 TEL: + 886 2 29993869 FAX: + 886 2 29993647
 Email: monitex@monitex.com.tw
www.monitex.com.tw

